





## Forceps & Hemostats

|           |       |         |
|-----------|-------|---------|
| Dumont    | Pages | 63 – 72 |
| Fine      | Pages | 73 – 80 |
| Moria     | Page  | 74      |
| S&T       | Pages | 75, 77  |
| Standard  | Pages | 81 – 91 |
| Hemostats | Pages | 92 – 97 |

*See page 189 for alloy and additional information.*



Dumont #5SF

The finest tips ever produced by Dumont.



Super fine tip  
0.025 x 0.005 mm  
Inox  
11 cm

No. 11252-00

Dumont #5  
Fine Tips

Finished under magnification and polished on the inside, the tips close flat to 3 mm and are less likely to open under excess pressure.



Biology tip  
0.05 x 0.01 mm  
Inox  
11 cm

No. 11254-20

Dumont #5  
Ceramic Coated

These forceps are non-reflective. The scratch and corrosion resistant coating guarantees durability.



Biology tip  
0.05 x 0.01 mm  
Ceramic coated inox  
11 cm

No. 11252-50

Dumont #5CO

These forceps have been honed down to make this the lightest and most precise Dumont forceps available.



Biology tip  
0.05 x 0.01 mm  
Dumostar  
11 cm

No. 11295-20

For information on  
Dumont alloys and tip styles,  
see page 189.

MOST  
POPULAR

Dumont Micro-Blunted, Atraumatic Tipped Forceps

The sides of these forceps are distinctively ribbed to provide a better holding surface for maximum control of the instrument. The shanks are arched toward each other to transmit more holding strength to the points. The smooth inside surfaces at the tips will close flat to 5 mm with gentle pressure, increasing the total holding area considerably. The microscopically matched tips are gently blunted to prevent tissue puncture and are recommended as micro needle holders.

FORCEPS



Straight  
0.10 x 0.06 mm  
Inox  
11 cm

Medical #5  
No. 11253-20




Off-center  
0.13 x 0.08 mm  
Inox  
11.5 cm

Medical #5A  
No. 11253-21



Short curve  
0.17 x 0.10 mm  
Inox  
12 cm

Medical #7S  
No. 11273-22



Curved  
0.17 x 0.10 mm  
Inox  
11.5 cm

Medical #7  
No. 11273-20



Angled 45°  
0.10 x 0.06 mm  
Inox  
11 cm

Medical #5/45  
No. 11253-25



Angled 90°  
0.10 x 0.06 mm  
Inox  
11 cm

Medical #5/90  
No. 11253-29

Dumont #55

Bevelled, like #5 but with extra light and fine shanks.



Biology tip  
0.05 x 0.02 mm  
11 cm

Inox  
No. 11255-20

Dumostar  
No. 11295-51

Dumont #5  
Mirror Finish Forceps

The last 5 mm on the inside of the tips are highly polished to help prevent tissue from sticking. Ideal for embryos, nerves, or any other adhesive tissue.



Standard tip  
0.10 x 0.06 mm  
Inox  
11 cm

No. 11251-23

Biology tip  
0.05 x 0.02 mm  
Inox  
11 cm

No. 11252-23

For information on  
Dumont alloys and tip  
styles, see page 189.

Dumont #5

These are the most widely used Dumont forceps.  
Bevelled, very fine shanks for best visibility.



Standard tip  
0.10 x 0.06 mm  
11 cm

Inox  
No. 11251-20

Dumoxel  
No. 11251-30

Dumostar  
No. 11295-00

Carbon  
0.08 x 0.04 mm  
**Not autoclavable**  
No. 11251-10

Biology tip  
0.05 x 0.02 mm  
11 cm

Inox  
No. 11252-20

Dumoxel  
No. 11252-30

Dumostar  
No. 11295-10

**TITANIUM**  
No. 11252-40

Dumont WA #5

These newly designed forceps offer improved stability with an ergonomic body that allows the forceps to fit more comfortably in your hand. A domed end with no soldering grooves allows the forceps to be thoroughly and easily cleaned, preventing biofilm residues.



Side View

Standard tip  
0.10 x 0.06 mm  
Inox  
12 cm

No. 11251-21

Biology tip  
0.05 x 0.02 mm  
Inox  
12 cm

No. 11252-21

For information on  
Dumont alloys and tip  
styles, see page 189.

Dumont #5XL



Standard tip  
0.10 x 0.06 mm  
Inox  
15 cm

No. 11253-10

**Dumont #L5**

Slide sleeve keeps tips together so specimens may be securely held and retained. Slide may be locked in place so that forceps can be used in regular fashion. **Autoclavable.**



Standard tip  
0.10 x 0.06 mm  
Inox  
11 cm

No. 11253-27

**Dumont #5  
Teflon® Coated**

These antimagnetic forceps are ideal for use in cryogenic environments and for EM samples.



Standard tip  
0.10 x 0.06 mm  
Teflon coated inox  
11 cm

No. 11626-11

**Dumont #5/45**

45° tips for improved visibility and control.



**MOST  
POPULAR**

Standard tip  
0.10 x 0.06 mm  
Dumoxel  
11 cm

No. 11251-35

**Dumont #5/45  
Cover Slip Forceps**

Designed to pick up cover slips out of multi-well plates. The tip width is squared for a better grip.



400%

1 mm from tip to bend  
Upper tip: 0.10 x 0.08 mm  
Lower tip: 0.10 x 0.02 mm  
Dumoxel  
11.5 cm

No. 11251-33

### Dumont #3

Tapered and sturdy shanks.



Standard tip  
0.17 x 0.10 mm  
12 cm

Dumoxel  
No. 11231-30

Dumostar  
No. 11293-00

### Dumont #3c



Standard tip  
0.17 x 0.10 mm  
Inox  
11 cm

No. 11231-20

### Dumont #6

Sharply angled tips for access to difficult areas.



Standard tip  
0.17 x 0.10 mm  
Inox  
11.5 cm

No. 11260-20

Dumont #7  
Fine Tips

Finished under magnification and polished on the inside, the tips close flat to 3 mm and therefore less likely to open under excess pressure.



Biology tip  
0.07 x 0.03 mm  
Inox  
11 cm  
No. 11274-20

For information on  
Dumont alloys and tip  
styles, see page 189.

Dumont #7



Standard tip  
0.17 x 0.10 mm  
11.5 cm  
Dumoxel  
No. 11271-30  
  
Dumostar  
No. 11297-00

•TITANIUM•

No. 11272-40

Dumont #7  
Ceramic Coated



Biology tip  
0.07 x 0.04 mm  
Ceramic coated inox  
11 cm  
No. 11272-50

Dumont #7B



Standard tip  
Serrated  
0.17 x 0.10 mm  
Inox  
11 cm  
No. 11270-20



**Dumont #4**

Similar to #5, but with broader shanks.



Standard tip  
0.13 x 0.08 mm  
Dumoxel  
11 cm

No. 11241-30

Biology tip  
0.06 x 0.04 mm  
Dumoxel  
11 cm

No. 11242-40

Dumostar  
No. 11294-00

**Dumont #2  
Laminectomy Forceps**

These forceps can be used to break off delicate bone during mouse laminectomies, as well as fine skull bone removal.



0.40 x 0.02 mm  
Inox  
12 cm

No. 11223-20

**Dumont #2AP  
Epoxy Coated**

Provides good insulation for up to 220 AC voltage work and highly impervious to chemical solutions. Frozen or super cold sections will not adhere to the gripping surfaces.



Tip width: 1.5 mm  
Epoxy coated Inox  
12 cm

No. 11220-21

**Dumont AA****Dumont AA  
Epoxy Coated**

Slip-proof acid and electric shock resistant epoxy coating covering 17.5 mm from the tips.



0.4 x 0.2 mm  
Inox  
12.5 cm

No. 11210-20



0.4 x 0.2 mm  
Epoxy coated inox  
12.5 cm

No. 11210-10

For information on  
Dumont alloys and tip  
styles, see page 189.

**Dumont SS**

The tips are smooth on the inside for a gentle grip. Applying full strength will not separate tips and a strong hold on the smallest specimen is ensured.



Standard tip  
0.20 x 0.12 mm  
Inox  
13.5 cm

No. 11203-23



Biology tip  
0.10 x 0.06 mm  
Dumoxel  
13.5 cm

No. 11200-33



Angled 45°

No. 11203-25

## Dumont Mini Forceps



|                                |                                |                                  |
|--------------------------------|--------------------------------|----------------------------------|
| 0.12 x 0.07 mm<br>Inox<br>7 cm | 0.14 x 0.06 mm<br>Inox<br>7 cm | 0.09 x 0.05 mm<br>Inox<br>8.2 cm |
| Style 3<br>No. 11200-10        | Style 4<br>No. 11200-12        | Style 5<br>No. 11200-14          |

## Fine Forceps Repair Kit

Realign and repair fine forceps with this convenient kit. Bent forceps tips can be restraightened with the chain nose pliers. The sharpening stone has two grooves, 1.6 mm and 0.8 mm, to assist in reshaping the forceps tips. O-rings clamp the forceps shut to ensure that both shanks are equally aligned and sharpened. Instrument oil assists in producing a smooth sharp edge. Detailed instructions are included.



Kit includes:

- 1 Sharpening stone for fine forceps with ten O-rings (No. 29008-22)
- 1 Instrument oil, 50 ml bottle (No. 29055-00)
- 1 Chain nose pliers

No. 29000-00

**MOST  
POPULAR**

## O-Rings for Forceps

These rubber O-rings fit our forceps and keep them closed with variable degrees of pressure on the tips.

Pkg. of 10  
No. 11200-00



**Fine Forceps Mirror Finish**

The last 5 mm on the inside of the tips are highly polished to help prevent tissue from sticking. Ideal for embryos, nerves, or other adhesive tissue.



0.10 x 0.06 mm  
11 cm

No. 11412-11      No. 11413-11

**Fine Forceps Broad Shanks**

Straight forceps with sturdy coarse tips.



0.2 x 0.12 mm  
12 cm

No. 11445-12

**Fine Forceps Self-Closing**

Ideal for lifting and holding EM grids. Consistent pressure at tips. Correctly tensioned "cross-over" shanks come to a microscopically matched point.



0.08 x 0.04 mm  
11.5 cm

No. 11480-11



0.10 x 0.06 mm  
11 cm

No. 11485-11



0.20 x 0.12 mm  
11.5 cm

No. 11487-11

These forceps have the most finely executed points available. They are extremely light, superbly balanced, and have a gentle spring action. The shanks are specially designed with an arch which ensures maximum force to the tips. The tips close flat for 6 mm from the points to give a large holding area. The tips have been expertly finished under magnification.



Moria



Moria



Moria

Tip width: 0.5 mm  
Blunt  
10 cm



Straight  
11 cm

Moria MC40  
No. 11370-40



Curved  
11 cm

Moria MC40/B  
No. 11370-42



Straight  
13 cm

Moria 9980  
No. 11399-80



Curved  
13 cm

Moria 9987  
No. 11399-87



Smooth  
Moria MC32/B  
No. 11373-22



Smooth  
Moria MC31/B  
No. 11373-12



Microerrated  
Moria MC32  
No. 11370-32



Microerrated  
Moria MC31  
No. 11370-31

### Delicate Moria Forceps

Balanced with gentle spring action and tips that are perfectly matched.


### SuperGrip Forceps

The tips are processed with MicroCrystals<sup>®</sup>, a method developed by S&T<sup>®</sup> to infuse tungsten carbide on stainless steel. The formation of this ultra-fine coating guarantees atraumatic handling of tissue while improving the grip dramatically. The recommended use is where forceps are utilized as needle holders or when slippery tissues might be encountered.



S&T<sup>®</sup>

0.3 x 0.25 mm  
11 cm



Straight

No. 00108-11



Angled 45°

No. 00109-11



S&T<sup>®</sup>

0.3 x 0.25 mm  
11 cm



The gripping surface is similar to a match head.



Straight

No. 00632-11




Angled 45°


No. 00649-11

**Bonn**

The 1x2 forceps are the finest of their kind available.




20x magnification




Smooth  
Tip width: 0.3 mm  
Tip platform: 5 mm  
7 cm

No. 11083-07




1x2 teeth  
Tip width: 0.45 mm  
Teeth width: 0.1 mm  
7 cm

No. 11084-07



20x magnification




Angled 45°  
Tip width: 0.4 mm  
Platform  
9 cm


No. 11063-07

**Mirror Finish  
Suture Tying Forceps**

A unique tip style fusing a suture tying platform and a 1x2 teeth tip for tissue manipulation.



20x magnification




Straight  
Tip width: 0.3 mm  
Platform with 1x2 teeth  
9 cm

No. 11090-09

**Suture Tying Forceps**

The platform tips close flat. This ensures a reliable grip on sutures.




Tip diameter: 0.3 mm  
Tying platform: 6 mm  
13.5 cm

No. 00272-13

**Suture Tying Forceps**

The platform tips close flat. This ensures a reliable grip on sutures.




Tip diameter: 0.5 mm  
Tying platform: 5 mm  
10 cm

No. 18025-10

**Round Handled Suture Tying Forceps**

Round handles allow easy fingertip adjustments.




Tip diameter: 0.35 mm  
Tying platform: 5.5 mm  
10.5 cm

No. 18026-10




**Vial Forceps**

Tips taper for 6 cm allowing easy access to specimens in vials.



Tip diameter: 0.5 mm  
15 cm

No. 11006-15

**Vetro Grasping Forceps**

Serrated  
Tip diameter: 0.7 mm  
14 cm

No. 11900-14



Tip width: 0.6 mm  
7 cm




Straight  
Serrated

No. 11064-07



Curved  
Serrated

No. 11065-07



Straight  
1x2 teeth

No. 11066-07



Curved  
1x2 teeth

No. 11067-07

## Bishop-Harman

Extra light, fenestrated handles ensure a good grip.




Tip width: 0.5 mm  
8.5 cm




Serrated

No. 11069-08



1x2 teeth

No. 11070-08



10x magnification

## Octagon Forceps

The octagon shaped grip provides a particularly secure finger hold.



Tip width: 0.8 mm  
9.5 cm




Straight  
Serrated

No. 11040-08



Curved  
Serrated

No. 11041-08



Straight  
1x2 teeth

No. 11042-08



Curved  
1x2 teeth

No. 11043-08

## Extra Fine Graefe Forceps



Tip width: 0.5 mm  
10 cm



Serrated  
No. 11150-10



Serrated  
No. 11151-10



Serrated  
No. 11152-10




Smooth  
No. 11158-10



Smooth  
No. 11159-10



Smooth  
No. 11160-10



1x2 teeth  
No. 11153-10



1x2 teeth  
No. 11154-10



1x2 teeth  
No. 11155-10

## Graefe Forceps



Tip width: See below  
10 cm



Serrated  
Tip width: 0.8 mm  
No. 11050-10



Serrated  
Tip width: 0.8 mm  
No. 11051-10

**MOST  
POPULAR**



Serrated  
Tip width: 0.8 mm  
No. 11052-10



Serrated  
Tip width: 0.8 mm  
No. 11049-10

**·TITANIUM·**

Tip width: 1.2 mm  
No. 11650-10


**·TITANIUM·**

Tip width: 1.2 mm  
No. 11651-10

**·TITANIUM·**

Tip width: 1.2 mm  
No. 11652-10

For more information  
see page 183.



1x2 teeth  
Tip width: 0.8 mm  
No. 11053-10



1x2 teeth  
Tip width: 0.8 mm  
No. 11054-10



1x2 teeth  
Tip width: 0.8 mm  
No. 11055-10

The tips of these forceps arch toward each other, closing flat. Exceptional gripping power is possible with minimal damage to tissue. A limit pin prevents excessive applied pressure.



Slide lock  
Smooth  
Tip width: 1 mm  
11 cm

No. 11062-10




Tip thickness: 0.2 mm  
10 cm

No. 11071-10

### Pierce Fixation Forceps

For exceptional gripping power with no damage to sclera, vascular sheath and epineurium.



Tissue retention  
hole of 0.3 mm  
behind tips.

Straight  
0.1 x 0.4 mm  
13 cm

No. 18155-13

**Plastic Forceps**

These forceps are non-magnetic, acid resistant and able to withstand temperatures up to 240°C (464°F). Tips do not spread under pressure.



Acetyl  
11.5 cm  
Pkg. of 10 forceps



Sharp straight  
0.40 x 0.25 mm

No. 11700-00



Blunt straight  
1.5 x 0.40 mm

No. 11700-03




Sharp curved  
0.30 x 0.20 mm

No. 11700-01



Blunt curved  
1.5 x 0.40 mm

No. 11700-04




Sharp 45°  
0.30 x 0.20 mm

No. 11700-02

**Cross Action Forceps**

These tissue transport forceps have a gentle action.




Cross serrated  
Tip width: 4 mm  
Tip depth: 8 mm  
10 cm

No. 11075-00

**Cover Glass Forceps**

10.5 cm



Straight

No. 11074-02

Angled

No. 11073-10

**Ring Forceps**

For transport of tissue and organs. With limit screw to prevent excessive applied pressure.



ID: 1 mm  
OD: 1.5 mm  
8.5 cm

No. 11101-09

ID: 2.2 mm  
OD: 3 mm  
9 cm

No. 11103-09

ID: 4.8 mm  
OD: 6 mm  
9 cm

No. 11106-09

**Circle Forceps**

Cupped forceps for atraumatic gripping of organs and tissues.



Tip diameter: 0.7 mm  
11.5 cm

No. 11072-11

**Live Insect Handling Forceps**

Very delicate.



Narrow  
10 cm

No. 26029-10

Wide  
10 cm


No. 26030-10

**Micro-Adson**

Our most delicate version of the popular Adson forceps.



Tip width: 0.8 mm  
12 cm



Serrated

1x2 teeth

No. 11018-12


No. 11019-12

**Micro-Adson with Fenestrated Handle**

Our delicate Micro-Adson forceps, with five windows to aid in stability and easy finger adjustments.



12 cm



Serrated

1x2 teeth

1x2 teeth with platform

Tip width: 0.8 mm  
No. 11020-12

Tip width: 0.5 mm  
No. 11032-12


Tip width: 0.8 mm  
No. 11022-12

**Adson**

Wide, strong grip with narrow tips.



Tip width: 1.5 mm  
12 cm



Serrated

1x2 teeth

Cross serrated tungsten carbide inlaid tip



No. 11006-12

No. 11027-12

No. 11506-12

**TUNGSTEN**  
**CARBIDE**

For more information  
see page 183.

**Adson DeBakey****Atraumatic  
DeBakey Serrations**


This tooth pattern, as shown, consists of one row of atraumatic teeth meshed with two opposite rows of atraumatic teeth.

These two forceps are ideal for gentle, atraumatic gripping of vessels, tissue and other delicate structures.



DeBakey serrations  
Tip width: 1.5 mm  
12.5 cm

No. 11617-12




DeBakey serrations  
Tip width: 1.5 mm  
15 cm

No. 11616-15

**Adson-Graefe**

Fine sharp teeth  
Tip width: 4.5 mm  
12.5 cm

No. 11030-12

**Adson-Brown**

Shark teeth  
Tip width: 1.4 mm  
12 cm

No. 11627-12



Slim pattern with long, narrow shanks and gentle spring action.



•TITANIUM•



•TITANIUM•

For more information  
see page 183.



Lighter spring action than standard pattern for grasping more delicate tissue.




•TITANIUM•


A pair of dental forceps, labeled F.S.T 11002-13, shown against a black background. The forceps have a long, slender handle with a textured grip and a curved, serrated tip. The label is printed in white at the bottom of the handle.

**Standard Pattern**

Broad pattern with strong spring action.



| Tip width (mm) | Length (cm) | Item         | Price |
|----------------|-------------|--------------|-------|
| 2.5            | 12          | No. 11000-12 |       |
| 2.8            | 13          | No. 11000-13 |       |
| 3.25           | 14.5        | No. 11000-14 |       |
| 3.3            | 16          | No. 11000-16 |       |
| 3.7            | 18          | No. 11000-18 |       |
| 3.75           | 20          | No. 11000-20 |       |
| 3.9            | 25          | No. 11000-25 |       |
| 4              | 30          | No. 11000-30 |       |



| Tip width (mm) | Length (cm) | Item         | Price |
|----------------|-------------|--------------|-------|
| 2.5            | 12          | No. 11001-12 |       |
| 2.8            | 13          | No. 11001-13 |       |
| 3.25           | 14.5        | No. 11001-14 |       |
| 3.3            | 16          | No. 11001-16 |       |
| 3.7            | 18          | No. 11001-18 |       |
| 3.75           | 20          | No. 11001-20 |       |

**Mouse Handling Forceps**

Long handles make picking up small rodents easy and help prevent bites.



Forceps with serrated silicone tips  
Gripping surface: 15 x 20 mm  
21 cm

No. 11036-20



Replacement silicone tips

Pkg of 2 tips  
No. 11036-90



Slim 1x2 teeth

10 cm  
No. 11023-10

14.5 cm  
No. 11023-14

12 cm  
No. 11023-12

15.5 cm  
No. 11023-15



1x2 teeth

12 cm  
No. 11021-12

14.5 cm  
No. 11021-14

15.5 cm  
No. 11021-15



2x3 teeth

14.5 cm  
No. 11022-14


15.5 cm  
No. 11022-15



1x2 teeth

14 cm  
No. 11029-14

Used with rongeur No. 16021-14  
(page 36) in laminectomy to securely  
hold spinal column.


**Gillies**

Serrated  
1x2 teeth  
15 cm

No. 11028-15

**Russian**

Special cupped pattern  
for superior grip.




15 cm

No. 11026-15

**Organ Holding Forceps**


With limit screw to prevent  
excessive applied pressure.



Serrated  
14.5 cm

No. 11031-15

## London Forceps



Angled  
Serrated  
16 cm

No. 11080-02

## Titanium Forceps



·TITANIUM·

For more information  
see page 183.



Straight  
Serrated  
16 cm


No. 11602-16

## Gerald Forceps



Straight  
Serrated  
18 cm

No. 11014-18




Curved  
Serrated  
18 cm

No. 11015-18

**Taylor Forceps**

Serrated and smooth versions available.

**Potts-Smith Forceps**

## Ultra Fine Hemostat



Tip width: 0.6 mm  
Smooth  
12 cm



Straight

No. 13020-12



Curved

No. 13021-12

## Fine Hemostat



Tip width: 1 mm  
Smooth  
12 cm



Straight

No. 13006-12



Curved

No. 13007-12

Micro-Mosquito



Tip width: 1.3 mm  
Serrated  
12 cm




Straight

No. 13010-12



Curved

No. 13011-12



Angled 90°

No. 13012-12

Halsted-Mosquito



Tip width: 1.1 mm  
Serrated  
12.5 cm



Straight

No. 13008-12



Curved

No. 13009-12



Straight

No. 13108-12



Curved

No. 13109-12





Crile




Tip width: 1.5 mm  
Serrated  
14 cm



Straight

No. 13004-14



Curved

No. 13005-14

Kelly




Tip width: 1.5 mm  
Serrated  
14 cm



Straight

No. 13018-14




Curved

No. 13019-14

Baby-Mixer



Tip width: 1.3 mm  
Serrated  
14 cm



Curved

No. 13013-14

Providence-Hospital



Tip width: 1.5 mm  
Serrated  
14 cm



Straight

No. 13014-14



Curved

No. 13015-14

Hartman



Tip width: 1 mm  
Serrated  
10 cm



Straight

No. 13002-10



Curved

No. 13003-10

MOST  
POPULAR

Backhaus Towel Forceps



9 cm

No. 11095-09



13 cm

No. 11095-13