



Scissors, Bone Instruments, Knives & Scalpels

Spring Scissors	Pages	3 – 14
Fine Scissors	Pages	15 – 28
Surgical Scissors	Pages	29 – 35
Bone Instruments	Pages	36 – 51
Knives & Scalpels	Pages	52 – 63



Micro Vannas

These Vannas-style spring scissors have the smallest blades ever produced. They are ideal where space is restricted and where absolute precision is demanded.



Straight
Sharp/Sharp

8 cm
Cutting edge: 2 mm
Tip Ø: 0,05 mm

No. 15000-03



Curved
Sharp/Sharp

8 cm
Cutting edge: 2 mm
Tip Ø: 0,05 mm

No. 15000-04

Micro Vannas

The short taper is ideal for small, precise cuts.



Straight
Sharp/Sharp

8 cm
Cutting edge: 2,5 mm
Tip Ø: 0,05 mm

No. 15000-08



Curved
Sharp/Sharp

8 cm
Cutting edge: 2,5 mm
Tip Ø: 0,05 mm

No. 15001-08



Angled to side
Sharp/Sharp

8 cm
Cutting edge: 2,5 mm
Tip Ø: 0,05 mm

No. 15002-08

Micro Vannas

The long taper is ideal for cutting vessel walls or when a longer cutting edge is needed.




15000-00



**FST
HAND
CRAFTED**

For more information
see page 38.



300%



Straight
Sharp/Sharp

8 cm
Cutting edge: 3 mm
Tip Ø: 0,05 mm

No. 15000-00



Curved
Sharp/Sharp

8 cm
Cutting edge: 3 mm
Tip Ø: 0,05 mm

No. 15000-10



**FST
HAND
CRAFTED**

For more information
see page 38.



Straight
Sharp/Sharp

10,5 cm
Cutting edge: 4 mm
Tip Ø: 0,15 mm

No. 15018-10



Curved
Sharp/Sharp

10,5 cm
Cutting edge: 4 mm
Tip Ø: 0,15 mm

No. 15019-10

Vannas-Tübingen



15003-08



**FST
HAND
CRAFTED**

For more information
see page 38.



Straight
Sharp/Sharp

8,5 cm
Cutting edge: 5 mm
Tip Ø: 0,075 mm

No. 15003-08



Angled to side
Sharp/Sharp

8,5 cm
Cutting edge: 5 mm
Tip Ø: 0,075 mm

No. 15008-08



Curved
Sharp/Sharp

8,5 cm
Cutting edge: 5 mm
Tip Ø: 0,075 mm

No. 15004-08



Angled up
Sharp/Sharp

9,5 cm
Cutting edge: 5 mm
Tip Ø: 0,075 mm

No. 15005-08

Vannas Scissors with Microserrations

Microserrated to minimize slippage.



FST
HAND
CRAFTED

For more information
see page 38.



Straight
Sharp/Sharp

8,5 cm
Cutting edge: 4 mm
Tip Ø: 0,1 mm

No. 15070-08



Straight
Sharp/Sharp

9 cm
Cutting edge: 5 mm
Tip Ø: 0,1 mm

No. 15007-08



Curved
Sharp/Sharp

8,5 cm
Cutting edge: 4 mm
Tip Ø: 0,1 mm

No. 15071-08

Miniature Vannas

Hand finished under high magnification, these tips are especially useful when incising vessel walls for catheter insertion.



Moria 9600



Straight
Sharp/Sharp

8 cm
Cutting edge: 5 mm
Tip Ø: 0,1 mm

No. 15396-00

Moria 9601



Curved
Sharp/Sharp

8 cm
Cutting edge: 5 mm
Tip Ø: 0,1 mm

No. 15396-01

Titanium Vannas

TITANIUM

For more information
see page 205.



**FST
HAND
CRAFTED**

For more information
see page 38.



Straight
Sharp/Sharp

9,5 cm
Cutting edge: 5 mm
Tip Ø: 0,3 mm

No. 15610-08



Curved
Sharp/Sharp

9,5 cm
5 mm
Tip Ø: 0,3 mm

No. 15615-08



Angled up
Sharp/Sharp

9,5 cm
5 mm
Tip Ø: 0,3 mm

No. 15617-08

Cohan-Vannas

Long cutting blades for
a more controlled cut.



**FST
HAND
CRAFTED**

For more information
see page 38.



Straight
Sharp/Sharp

6 cm
Cutting edge: 6 mm
Tip Ø: 0,1 mm

No. 15000-02



Curved
Sharp/Sharp

6 cm
6 mm
Tip Ø: 0,1 mm

No. 15000-01

Straight
Blunt/Blunt

6 cm
6 mm
Tip Ø: 0,1 mm

No. 15000-12

Curved
Blunt/Blunt

6 cm
6 mm
Tip Ø: 0,1 mm

No. 15000-11

ToughCut® Spring Scissors

**TOUGH
CUT**

For more information
see page 205.



**FST
HAND
CRAFTED**

For more information
see page 38.



Straight
Sharp/Sharp

12,5 cm
6 mm
Tip Ø: 0,3 mm

No. 15124-12



Curved
Sharp/Sharp

12,5 cm
6 mm
Tip Ø: 0,3 mm

No. 15123-12

Round Handled Vannas

Round handles allow easy fingertip adjustments.



For more information
see page 38.



Straight
Sharp/Sharp

12,5 cm
Cutting edge: 6 mm
Tip Ø: 0,075 mm

No. 15400-12




Curved
Sharp/Sharp

12,5 cm
Cutting edge: 6 mm
Tip Ø: 0,075 mm

No. 15401-12



For more information
see page 38.



Straight
Sharp/Sharp

14 cm
Cutting edge: 6 mm
Tip Ø: 0,15 mm

No. 15020-15



Curved
Sharp/Sharp

14 cm
Cutting edge: 6 mm
Tip Ø: 0,15 mm

No. 15021-15



For more information about
our ceramic instruments
please visit page 203.



For more information
see page 38.



Straight
Sharp/Sharp

10 cm
Cutting edge: 7 mm
Tip Ø: 0,1 mm

No. 15750-11



Curved
Sharp/Sharp

10 cm
Cutting edge: 7 mm
Tip Ø: 0,1 mm

No. 15751-11

Vannas-Wolff



Moria MC52



Straight
Sharp/Sharp


10 cm

Cutting edge: 7 mm

Tip Ø: 0,15 mm

No. 15370-52

Moria MC50



Curved
Sharp/Sharp

10 cm

Cutting edge: 7 mm

Tip Ø: 0,15 mm

No. 15370-50

Vannas

Ultra long cutting edge
without sacrificing tip width.



**FST
HAND
CRAFTED**

For more information
see page 38.



Straight
Sharp/Sharp

8,5 cm

Cutting edge: 8 mm

Tip Ø: 0,1 mm

No. 15009-08

Castroviejo



**FST
HAND
CRAFTED**

For more information
see page 38.



Angled to side
Sharp/Sharp

10 cm

Cutting edge: 10 mm

Tip Ø: 0,125 mm

No. 15006-09

Castroviejo with Ball Tip

Same as No. 15006-09 but with a ball tip to prevent piercing of underlying tissue.



**FST
HAND
CRAFTED**

For more information
see page 38.



Angled to side
Sharp/Sharp

9 cm
Cutting edge: 8 mm
Tip Ø: 1,3 mm Ball

No. 15033-09

Spring Scissors



**FST
HAND
CRAFTED**

For more information
see page 38.



Straight
Sharp/Sharp

10 cm
Cutting edge: 8 mm
Tip Ø: 0,15 mm

No. 15024-10



Straight
Blunt/Blunt

10 cm
Cutting edge: 8 mm
Tip Ø: 0,15 mm

No. 15025-10



Curved
Sharp/Sharp

10 cm
Cutting edge: 8 mm
Tip Ø: 0,15 mm

No. 15023-10



**FST
HAND
CRAFTED**


For more information
see page 38.



Straight
Sharp/Sharp

11,5 cm
Cutting edge: 10 mm
Tip Ø: 0,1 mm

No. 15040-11



Angled up
Sharp/Sharp

11,5 cm
Cutting edge: 10 mm
Tip Ø: 0,1 mm

No. 15042-11

Noyes



For more information
see page 38.



Straight
Sharp/Sharp

12 cm
Cutting edge: 14 mm
Tip Ø: 0,2 mm


No. 15012-12



Curved
Sharp/Sharp

12 cm
Cutting edge: 14 mm
Tip Ø: 0,2 mm

No. 15011-12



Angled up
Sharp/Sharp

12 cm
Cutting edge: 14 mm
Tip Ø: 0,2 mm

No. 15013-12

Tungsten Carbide Noyes

Tungsten carbide cutting edges
provide long lasting sharpness.



For more information
see page 203.



For more information
see page 38.



Straight
Sharp/Sharp

12 cm
Cutting edge: 14 mm
Tip Ø: 0,275 mm

No. 15514-12

Castroviejo



For more information
see page 38.



Sharply curved
Sharp/Sharp

10 cm
Cutting edge: 10 mm
Tip Ø: 0,2 mm

No. 15017-10

Ceramic Coated Titanium Castroviejo

Westcott

Fine Angled Spring Scissors

These spring scissors aid in cutting tissues from various angles.



For more information
see page 38.



For more information
see page 203.



For more information about
our ceramic instruments
please visit page 203.



Straight
Sharp/Sharp

11 cm
Cutting edge: 5 mm
Tip Ø: 0,2 mm

No. 15650-11



Curved
Sharp/Sharp

11 cm
Cutting edge: 5 mm
Tip Ø: 0,2 mm

No. 15651-11



For more information
see page 38.



Slightly curved
Sharp/Sharp

11 cm
Cutting edge: 15 mm
Tip Ø: 0,25 mm

No. 15015-11



For more information
see page 38.



Angled up
Sharp/Sharp

9,5 cm
Cutting edge: 2 mm
Tip Ø: 0,1 mm

No. 15010-09



Angled up
Sharp/Sharp

10 cm
Cutting edge: 4 mm
Tip Ø: 0,125 mm

No. 15010-10

Trident II

A re-release of a popular spring scissors now with a new innovative design for a better feel.



Straight
Sharp/Sharp

8 cm
Cutting edge: 8 mm
Tip Ø: 0,2 mm

No. 15110-08



Angled to side
Sharp/Sharp

8 cm
Cutting edge: 8 mm
Tip Ø: 0,2 mm


No. 15111-08

Pascheff-Wolff

Very delicate blades.



Moria MC19/B



Angled up
Sharp/Sharp

10,5 cm
Cutting edge: 3 mm
Tip Ø: 0,15 mm

No. 15371-92

Long Handle Spring Scissors



Straight
Sharp/Sharp

18,5 cm
Cutting edge: 5 mm
Tip Ø: 0,2 mm

No. 15028-18



Curved
Sharp/Sharp

18,5 cm
Cutting edge: 5 mm
Tip Ø: 0,2 mm

No. 15027-18


Micro Scissors

With microserrations to minimize slippage.




**FST
HAND
CRAFTED**

For more information
see page 38.



300%



300%



Straight
Sharp/Sharp

14,5 cm
Cutting edge: 2,5 mm
Tip Ø: 0,175 mm

No. 15030-14



Angled to side
Sharp/Sharp

14,5 cm
Cutting edge: 2,5 mm
Tip Ø: 0,175 mm

No. 15031-14


Vitrectomy Scissors

The finest blades on spring scissors. Designed for ophthalmic surgery (vitrectomy). The spring handles are made of ribbed plastic to ensure a solid grip. **Autoclavable.**



**FST
HAND
CRAFTED**


For more information
see page 38.



300%



300%



Straight
Sharp/Sharp

14 cm
Cutting edge: 2,2 mm
Tip Ø: 0,15 mm

No. 15036-14



Angled 90°
Sharp/Sharp

14 cm
Cutting edge: 1 mm
Tip Ø: 0,1 mm

No. 15037-14

Ultra Fine Clipper Scissors II

For microscopic dissections of soft and mucous tissue. The minute fineness of the blades allows clear visualization of the tissue being dissected.



Uses:

- Removing sheath of ganglia to access neurons in aplysia.
- Cutting single nerve axons in aplysia, mollusk and snail.
- Stripping fascia off arteries and muscles.
- Ultra fine brain tissue dissection.
- Cutting of minute fiber optic strands.

Straight
Sharp/Sharp

10 cm
Cutting edge: 1,5–2 mm
Tip Ø: 0,05–0,08 mm
Blade shape: Triangular
Blade thickness: 0,1 mm

No. 15300-00

Biemer Scissors

These scissors are designed to securely cut vessels without slipping.



Straight
Blunt/Blunt

8,5 cm
Cutting edge: 4 mm
Tip Ø: 2 mm

No. 15032-08

Straight
Blunt/Blunt

13 cm
Cutting edge: 2 mm
Tip Ø: 2 mm

No. 15032-13

Extra Fine Bonn

Our finest iris scissors, best suited for microdissection under high magnification.

**MOST
POPULAR**

14084-08
14085-08



Straight
Sharp/Sharp

8,5 cm
Cutting edge: 13 mm

No. 14084-08



Curved
Sharp/Sharp

8,5 cm
Cutting edge: 13 mm

No. 14085-08



Straight
Sharp/Blunt

8,5 cm
Cutting edge: 13 mm

No. 14083-08

Bonn

Highly polished.



Straight
Sharp/Sharp

9 cm
Cutting edge: 13 mm

No. 14184-09



Curved
Sharp/Sharp

9 cm
Cutting edge: 13 mm

No. 14185-09

ToughCut® Bonn-Strabismus

A fusion of two light and delicate iris scissors with ToughCut cutting edges.

**TOUGH
CUT**

For more information
see page 203.



Straight
Blunt/Blunt

9 cm
Cutting edge: 14 mm

No. 14084-09



Curved
Blunt/Blunt

9 cm
Cutting edge: 14 mm

No. 14085-09

Hardened Fine Iris Scissors

These scissors have an anti-reflective finish to reduce glare.

**MOST
POPULAR**

14090-09



Straight
Sharp/Sharp

9 cm
Cutting edge: 24 mm

No. 14090-09



Curved
Sharp/Sharp

9 cm
Cutting edge: 24 mm

No. 14091-09

11 cm
Cutting edge: 26 mm

No. 14090-11

11 cm
Cutting edge: 26 mm

No. 14091-11

Fine Iris Scissors

Made from Martensitic stainless steel, combined with Molybdenum and Vanadium.



Straight
Sharp/Sharp

10,5 cm
Cutting edge: 26 mm

No. 14094-11



Curved
Sharp/Sharp

10,5 cm
Cutting edge: 26 mm

No. 14095-11

Iris Scissors

Delicate pattern.

**MOST
POPULAR**


14060-11

Straight
Sharp/Sharp9 cm
Cutting edge: 21 mm


No. 14060-09

Curved
Sharp/Sharp9 cm
Cutting edge: 21 mm

No. 14061-09

Curved to side
Sharp/Sharp9 cm
Cutting edge: 21 mm

No. 14062-09

Angled to side
Sharp/Sharp9 cm
Cutting edge: 21 mm

No. 14063-09

10,5 cm
Cutting edge: 23 mm

No. 14060-10

10,5 cm
Cutting edge: 23 mm

No. 14061-10

10,5 cm
Cutting edge: 23 mm

No. 14062-10

11,5 cm
Cutting edge: 23 mm

No. 14063-11

11,5 cm
Cutting edge: 25 mm

No. 14060-11

11,5 cm
Cutting edge: 25 mm

No. 14061-11

Left Handed Iris Scissors**LEFT
HANDED** Straight
Sharp/Sharp11,5 cm
Cutting edge: 25 mm

No. 14160-10

Curved
Sharp/Sharp11,5 cm
Cutting edge: 25 mm

No. 14161-10

Moria Iris Scissors



Moria 4877



Straight
Sharp/Sharp

12,5 cm
Cutting edge: 25 mm

No. 14370-22

Moria 4878



Shallow curve
Sharp/Sharp

10,5 cm
Cutting edge: 25 mm

No. 14370-23

Extra Narrow Iris Scissors

Extra narrow, delicate pattern. Very
slim blades with gentle cutting action.




Side-view showing
extra narrow pattern

Straight
Sharp/Sharp

10,5 cm
Cutting edge: 23 mm

No. 14088-10

Large Loops Scissors



Straight
Sharp/Sharp

10 cm
Cutting edge: 25 mm

No. 14040-10



Curved
Sharp/Sharp

10 cm
Cutting edge: 25 mm

No. 14041-10

Wagner



14070-12
14071-12



Straight
Sharp/Sharp

12 cm
Cutting edge: 30 mm

No. 14068-12



Curved
Sharp/Sharp

12 cm
Cutting edge: 30 mm

No. 14069-12



Straight
Sharp/Sharp

12 cm
Cutting edge: 30 mm

No. 14070-12



Curved
Sharp/Sharp

12 cm
Cutting edge: 30 mm

No. 14071-12

ToughCut®

MOST
POPULAR14058-09
14058-11TOUGH
CUTFor more information
see page 203.Straight
Sharp/Sharp9 cm
Cutting edge: 22 mm

No. 14058-09

Curved
Sharp/Sharp9 cm
Cutting edge: 22 mm

No. 14059-09

11,5 cm
Cutting edge: 25 mm

No. 14058-11

11,5 cm
Cutting edge: 25 mm

No. 14059-11

Tungsten Carbide

MOST
POPULAR14568-09
14568-12TUNGSTEN
CARBIDEFor more information
see page 203.Straight
Sharp/Sharp9 cm
Cutting edge: 22 mm

No. 14568-09

Curved
Sharp/Sharp9 cm
Cutting edge: 22 mm

No. 14569-09

11,5 cm
Cutting edge: 25 mm

No. 14568-12

11,5 cm
Cutting edge: 25 mm

No. 14569-12

ToughCut® and Tungsten Carbide

For more information
see page 203.



For more information
see page 203.



Straight
Sharp/Sharp

9 cm
Cutting edge: 22 mm

No. 14558-09



Curved
Sharp/Sharp

9 cm
Cutting edge: 22 mm

No. 14559-09

11,5 cm
Cutting edge: 25 mm

No. 14558-11

11,5 cm
Cutting edge: 25 mm

No. 14559-11

CeramaCut®

These ceramic-coated scissors are the cleanest cutting, longest lasting scissors ever produced. Utilizing the unique properties of a ceramic coating, they will last six times longer than tungsten carbide instruments and cut even better than ToughCut® scissors. The black, corrosion-resistant coating has a non-glare finish.



For more information
see page 203.



Straight
Sharp/Sharp

9 cm
Cutting edge: 22 mm

No. 14958-09



Curved
Sharp/Sharp

9 cm
Cutting edge: 22 mm

No. 14959-09

11,5 cm
Cutting edge: 25 mm

No. 14958-11

11,5 cm
Cutting edge: 25 mm

No. 14959-11

Fine Scissors

Lexer-Baby



Straight
Sharp/Blunt

10 cm
Cutting edge: 27 mm

No. 14028-10



Curved
Sharp/Blunt

10 cm
Cutting edge: 27 mm

No. 14029-10



Straight
Blunt/Blunt

10 cm
Cutting edge: 25 mm

No. 14078-10



Curved
Blunt/Blunt

10 cm
Cutting edge: 25 mm

No. 14079-10

Strabismus

Tungsten Carbide Strabismus

**TUNGSTEN
CARBIDE**

For more information
see page 203.



Straight
Blunt/Blunt

9 cm
Cutting edge: 20 mm

No. 14074-09

11,5 cm
Cutting edge: 23 mm

No. 14074-11



Curved
Blunt/Blunt

9 cm
Cutting edge: 20 mm

No. 14075-09

11,5 cm
Cutting edge: 23 mm

No. 14075-11



Straight
Blunt/Blunt

9 cm
Cutting edge: 20 mm

No. 14574-09

11 cm
Cutting edge: 23 mm

No. 14574-11



Curved
Blunt/Blunt

9 cm
Cutting edge: 20 mm

No. 14575-09

11 cm
Cutting edge: 23 mm

No. 14575-11

Walton



Slightly curved
Sharp/Sharp

9,5 cm
Cutting edge: 15 mm

No. 14077-09



Slightly curved
Sharp/Sharp

10,5 cm
Cutting edge: 18 mm

No. 14077-10

Stevens Tenotomy



Straight
Sharp/Sharp

11 cm
Cutting edge: 21 mm

No. 14064-11



Straight
Blunt/Blunt

11 cm
Cutting edge: 21 mm

No. 14066-11



Curved
Sharp/Sharp

11 cm
Cutting edge: 21 mm

No. 14065-11



Curved
Blunt/Blunt

11 cm
Cutting edge: 21 mm

No. 14067-11

LaGrange

One edge is serrated to minimize slippage of tissue during cutting. Curved handle helps when cutting from an awkward angle.



Side View



Curved
Sharp/Sharp

11,5 cm
Cutting edge: 18 mm

No. 14173-12

Shea

Transplant scissors with light and delicate pattern.




Curved
Blunt/Blunt

12 cm
Cutting edge: 20 mm

No. 14105-12

Dissector Scissors

Slim blades for cutting tissue.



Angled to side
Sharp/Sharp

10 cm
Cutting edge: 8 mm

No. 14081-09

Dissector Scissors

Similar to 14081-09, but with stronger blades for cutting fine bone.



Angled to side
Sharp/Sharp

10,5 cm
Cutting edge: 12 mm

No. 14082-09

Spencer Ligature

Used for cutting stitches or suture.



Straight
Sharp/Sharp

9 cm
Cutting edge: 20 mm

No. 14076-09

Bonn Artery

Iris scissors with one ball tip to prevent piercing of underlying tissue.



Straight
Sharp/Ball tip

9 cm
Cutting edge: 13 mm
Tip Ø: 1,5 mm Ball

No. 14086-09



Straight
Sharp/Ball tip

11,5 cm
Cutting edge: 25 mm
Tip Ø: 1,5 mm Ball

No. 14080-11

Tissue Separating Scissors

Straight
Blunt/Blunt


11,5 cm
Cutting edge: 23 mm

No. 14072-10

Graefe**Adjustable Loops Fine Scissors****Kaye**

The finger loops can be enlarged to fit almost any size.

These scissors have semi-sharp outer edges to better facilitate dissection.



Sharply angled up
Sharp/Sharp

10,5 cm
Cutting edge: 14 mm

No. 14037-10



Straight
Sharp/Sharp

9,5 cm
Cutting edge: 14 mm

No. 14106-09



Straight
Blunt/Blunt

9,5 cm
Cutting edge: 14 mm

No. 14108-09



Slightly curved
Blunt/Blunt

10,5 cm
Cutting edge: 22 mm

No. 14171-10

Standard Pattern

MOST
POPULAR

14001-12 – 14001-18

LEFT
HANDED

14008-14

Our standard surgical scissors are recognized for these important features:

- Excellent feel and handling capability.
- Overall execution that marks them as fine examples of German craftsmanship.



Length (cm)	Cutting edge (mm)	Straight Sharp/Sharp	Straight Sharp/Blunt	Straight Blunt/Blunt
12	40	No. 14002-12	No. 14001-12	No. 14000-12
13	40	No. 14002-13	No. 14001-13	No. 14000-13
14,5	45	No. 14002-14	No. 14001-14	No. 14000-14
16,5	50	No. 14002-16	No. 14001-16	No. 14000-16
18,5	55		No. 14001-18	No. 14000-18
20	60			No. 14000-20

Length (cm)	Cutting edge (mm)	Straight Sharp/Blunt
14,5	45	No. 14008-14



Length (cm)	Cutting edge (mm)	Curved Sharp/Sharp	Curved Sharp/Blunt	Curved Blunt/Blunt
12	40	No. 14005-12	No. 14004-12	No. 14003-12
13	40	No. 14005-13	No. 14004-13	No. 14003-13
14,5	45	No. 14005-14	No. 14004-14	No. 14003-14
16,5	50	No. 14005-16	No. 14004-16	No. 14003-16
18,5	55		No. 14004-18	No. 14003-18
20	60			No. 14003-20

ToughCut® Standard Pattern



14054-13



For more information
see page 203.



80% of
actual size



Straight
Sharp/Blunt

13 cm
Cutting edge: 40 mm

No. 14054-13



Straight
Sharp/Blunt

17,5 cm
Cutting edge: 60 mm

No. 14130-17

Tungsten Carbide Standard Pattern



For more information
see page 203.



80% of
actual size



Straight
Sharp/Sharp

14,5 cm
45 mm

No. 14502-14



Straight
Sharp/Blunt

14,5 cm
45 mm

No. 14501-14



Curved
Sharp/Sharp

14,5 cm
45 mm

No. 14505-14



Curved
Sharp/Blunt

14,5 cm
45 mm

No. 14504-14

Standard Pattern with Serrations



80% of
actual size



Straight
Sharp/Blunt

14 cm
40 mm

No. 14007-14

**Standard Pattern with
Large Finger Loops**



80% of
actual size



Straight
Sharp/Blunt

14,5 cm
Cutting edge: 45 mm

No. 14101-14

**Tungsten Carbide
Mini-Lister**



For more information
see page 203.



Angled to side
Blunt/Blunt

9 cm
Cutting edge: 20 mm

No. 14531-09

Lister

Excellent skin cutting scissors. Longer, blunt
blade prevents damaging tissue underneath
skin without losing the cutting line.



80% of
actual size




Angled to side
Blunt/Blunt

14 cm
Cutting edge: 30 mm

No. 14131-14

Metzenbaum Baby



80% of
actual size



Straight
Blunt/Blunt

13 cm
Cutting edge: 30 mm

No. 14018-13



Curved
Blunt/Blunt

13 cm
30 mm

No. 14019-13

14,5 cm
Cutting edge: 35 mm

No. 14018-14

14,5 cm
35 mm

No. 14019-14

18 cm
Cutting edge: 40 mm

No. 14018-18

18 cm
40 mm

No. 14019-18

Tungsten Carbide Metzenbaum Fino

**TUNGSTEN
CARBIDE**

For more information
see page 203.



80% of
actual size



Straight
Blunt/Blunt
Ultra fine tips

14,5 cm
35 mm

No. 14508-15



Straight
Blunt/Blunt

14,5 cm
35 mm

No. 14518-14



Curved
Blunt/Blunt
Ultra fine tips

14,5 cm
35 mm

No. 14509-15



Curved
Blunt/Blunt

14,5 cm
35 mm

No. 14519-14

18 cm
40 mm

No. 14518-18

18 cm
40 mm

No. 14519-18


Metzenbaum

ToughCut® Metzenbaum



For more information
see page 203.

SCISSORS



80% of
actual size



80% of
actual size



Straight
Sharp/Blunt

14,5 cm
Cutting edge: 30 mm

No. 14024-14



Straight
Blunt/Blunt

14,5 cm
Cutting edge: 30 mm

No. 14016-14



Curved
Blunt/Blunt

14,5 cm
Cutting edge: 30 mm

No. 14017-14

18 cm
Cutting edge: 35 mm

No. 14016-18

18 cm
Cutting edge: 35 mm

No. 14017-18



Straight
Blunt/Blunt

14,5 cm
Cutting edge: 30 mm

No. 14116-14



Curved
Blunt/Blunt

14,5 cm
Cutting edge: 30 mm

No. 14117-14

18 cm
Cutting edge: 35 mm

No. 14116-18

18 cm
Cutting edge: 35 mm

No. 14117-18

Mayo

ToughCut® Mayo

Tungsten Carbide Mayo



For more information
see page 203.



For more information
see page 203.



14512-15 & 14512-17
For more information
see page 203.



80% of
actual size



80% of
actual size



80% of
actual size



Straight
Blunt/Blunt

15 cm
Cutting edge: 40 mm

No. 14010-15



Curved
Blunt/Blunt

15 cm
40 mm

No. 14011-15

17 cm
Cutting edge: 45 mm

No. 14010-17

17 cm
45 mm

No. 14011-17



Straight
Blunt/Blunt

15 cm
40 mm

No. 14110-15



Curved
Blunt/Blunt

15 cm
40 mm

No. 14111-15

17 cm
45 mm

No. 14110-17

17 cm
45 mm

No. 14111-17



Straight
Blunt/Blunt

15 cm
40 mm

No. 14510-15



Curved
Blunt/Blunt

15 cm
40 mm

No. 14511-15



Straight
Blunt/Blunt

15 cm
40 mm

No. 14512-15

17 cm
45 mm

No. 14510-17

17 cm
45 mm

No. 14511-17


17 cm
45 mm

No. 14512-17

Mayo-Stille



80% of
actual size



Straight
Blunt/Blunt

15 cm
Cutting edge: 40 mm

No. 14012-15



Curved
Blunt/Blunt

15 cm
Cutting edge: 40 mm

No. 14013-15

17 cm
Cutting edge: 45 mm

No. 14012-17

17 cm
Cutting edge: 45 mm

No. 14013-17

Mayo-Noble



80% of
actual size



Straight
Blunt/Blunt

16,5 cm
Cutting edge: 45 mm

No. 14014-17



Curved
Blunt/Blunt

16,5 cm
Cutting edge: 45 mm

No. 14015-17

Decapitation Scissors

Stop guard at tips to prevent slippage.



70% of
actual size



Straight
Blunt/Blunt

21 cm
Cutting edge: 55 mm

No. 14200-21