




Forceps & Hemostats

Dumont	Pages	65 – 74
Fine	Pages	75 – 84
Moria	Page	78
S&T	Pages	79, 81
Standard	Pages	85 – 95
Hemostats	Pages	96 – 101

See page 209 for alloy and additional information.



Dumont #5SF

The finest tips ever produced by Dumont.



Super fine tip

11 cm
Inox
0,025 x 0,005 mm

No. 11252-00

Dumont #5 Fine Tips

Finished under magnification and polished on the inside, the tips close flat to 3 mm and are less likely to open under excess pressure.



Biology tip

11 cm
Inox
0,05 x 0,01 mm

No. 11254-20

Dumont #5 Ceramic Coated

These forceps are non-reflective. The scratch and corrosion resistant coating guarantees durability.



For more information about our ceramic instruments please visit page 203.



Biology tip

11 cm
Ceramic coated inox
0,05 x 0,01 mm

No. 11252-50

Dumont #5C0

These forceps have been honed down to make this the lightest and most precise Dumont forceps available.



Biology tip

11 cm
Dumostar
0,05 x 0,01 mm

No. 11295-20

FORCEPS

For information on Dumont alloys and tip styles, see page 209.

Dumont Micro-Blunted, Atraumatic Tipped Forceps

The sides of these forceps are distinctively ribbed to provide a better holding surface for maximum control of the instrument. The shanks are arched toward each other to transmit more holding strength to the points. The smooth inside surfaces at the tips will close flat to 5 mm with gentle pressure, increasing the total holding area considerably. The microscopically matched tips are gently blunted to prevent tissue puncture and are recommended as micro needle holders.

FORCEPS

For information on
Dumont alloys
and tip styles,
see page 209.



Medical #5

Straight

11 cm
Inox
0,10 x 0,06 mm

No. 11253-20



Medical #5A

Off-center

11,5 cm
Inox
0,13 x 0,08 mm

No. 11253-21




Medical #7S

Short curve

12 cm
Inox
0,17 x 0,10 mm

No. 11273-22



Medical #7

Curved

11,5 cm
Inox
0,17 x 0,10 mm

No. 11273-20



Medical #5/45

Angled 45°

11 cm
Inox
0,10 x 0,06 mm

No. 11253-25



Medical #5/90

Angled 90°

11 cm
Inox
0,10 x 0,06 mm

No. 11253-29

Dumont #55

Bevelled, like #5 but with extra light and fine shanks.



Biology tip

11 cm
Inox
0,05 x 0,02 mm
No. 11255-20

11 cm
Dumostar
0,05 x 0,02 mm
No. 11295-51

Dumont #5 Mirror Finish Forceps

The last 5 mm on the inside of the tips are highly polished to help prevent tissue from sticking. Ideal for embryos, nerves, or any other adhesive tissue.



Standard tip

11 cm
Inox
0,10 x 0,06 mm
No. 11251-23



Biology tip

11 cm
Inox
0,05 x 0,02 mm
No. 11252-23

Dumont #5

These are the most widely used Dumont forceps. Bevelled, very fine shanks for best visibility.



11251-20
11252-20



11251-10



11252-40
For more information
see page 203.



Standard tip

11 cm
Inox
0,10 x 0,06 mm
No. 11251-20

Dumoxel

No. 11251-30

Dumostar

No. 11295-00

11 cm
Carbon
0,08 x 0,04 mm
No. 11251-10



Biology tip

11 cm
Inox
0,05 x 0,02 mm
No. 11252-20

Dumoxel

No. 11252-30

Dumostar

No. 11295-10

Titanium
No. 11252-40

For information on
Dumont alloys
and tip styles,
see page 209.

Dumont WA #5

These specially designed forceps offer improved stability with an ergonomic body that allows the forceps to fit more comfortably in your hand. A domed end with no soldering grooves allows the forceps to be thoroughly and easily cleaned, preventing biofilm residues.

For information on
Dumont alloys
and tip styles,
see page 209.



Side View



Standard tip

12 cm
Inox
0,10 x 0,06 mm

No. 11251-21



Biology tip

12 cm
Inox
0,05 x 0,02 mm

No. 11252-21

Dumont #5XL

Standard tip

15 cm
Inox
0,10 x 0,06 mm

No. 11253-10

Dumont #L5

Slide sleeve keeps tips together so specimens may be securely held and retained. Slide may be locked in place so that forceps can be used in regular fashion.



Standard tip
11 cm
Inox
0,10 x 0,06 mm

No. 11253-27

Dumont #5 Teflon® Coated

These antimagnetic forceps are ideal for use in cryogenic environments and for EM samples. Teflon® coating provides a non-stick surface.



Standard tip
11 cm
Teflon coated inox
0,10 x 0,06 mm

No. 11626-11

Dumont #5/45

45° tips for improved visibility and control.

MOST
POPULAR



Standard tip
11 cm
Dumoxel
0,10 x 0,06 mm

No. 11251-35

Dumont #5/45 Cover Slip Forceps

Designed to pick up cover slips out of multi-well plates. The tip width is squared for a better grip.



300%



1 mm from tip to bend
11,5 cm
Dumoxel
Upper tip: 0,10 x 0,08 mm
Lower tip: 0,10 x 0,02 mm

No. 11251-33

For information on
Dumont alloys
and tip styles,
see page 209.

FORCEPS

Dumont #3

Tapered and
sturdy shanks.



Standard tip

12 cm
Dumoxel
0,17 x 0,10 mm

No. 11231-30



Standard tip

12 cm
Dumostar
0,17 x 0,10 mm

No. 11293-00

Dumont #3c

Standard tip

11 cm
Inox
0,17 x 0,10 mm

No. 11231-20

Dumont #6

Sharply angled tips for
access to difficult areas.



Standard tip

11,5 cm
Inox
0,17 x 0,10 mm

No. 11260-20

Dumont #7 Fine Tips

Finished under magnification and polished on the inside, the tips close flat to 3 mm and therefore less likely to open under excess pressure.



Biology tip
11 cm
Inox
0,07 x 0,03 mm
No. 11274-20

Dumont #7

•TITANIUM•

11272-40
For more information
see page 203.



Standard tip
Dumoxel
No. 11271-30

Dumostar
No. 11297-00

Biology tip
11,5 cm
Dumoxel
0,07 x 0,04 mm
No. 11272-30

Dumostar
No. 11297-10

Titanium
No. 11272-40

Dumont #7 Ceramic Coated

CERAMIC
COATED

For more information about
our ceramic instruments
please visit page 203.



Biology tip
11 cm
Ceramic coated inox
0,07 x 0,04 mm
No. 11272-50

Dumont #7B

SERRATED



Standard tip
11 cm
Inox
0,17 x 0,10 mm
No. 11270-20

For information on
Dumont alloys
and tip styles,
see page 209.

FORCEPS

Dumont #4

Similar to #5, but
with broader shanks.



Standard tip

11 cm
Dumoxel
0,13 x 0,08 mm

No. 11241-30

Dumostar

No. 11294-00



Biology tip

11 cm
Dumoxel
0,06 x 0,04 mm

No. 11242-40

Dumont #2 Laminectomy Forceps

These forceps can be used to break off
delicate bone during mouse laminectomies,
as well as fine skull bone removal.




Standard tip

12 cm
Inox
0,40 x 0,02 mm

No. 11223-20

Dumont #2AP Epoxy Coated

Provides good insulation for up to 220 AC voltage
work and highly impervious to chemical
solutions. Frozen or super cold sections will
not adhere to the gripping surfaces.



Standard tip

12 cm
Epoxy coated Inox
Tip width: 1,5 mm

No. 11220-21

For information on
Dumont alloys
and tip styles,
see page 209.

Dumont AA

Dumont AA Epoxy Coated

Slip-proof acid and electric shock resistant epoxy coating covering 17,5 mm from the tips.

Dumont SS

The tips are smooth on the inside for a gentle grip. Applying full strength will not separate tips and a strong hold on the smallest specimen is ensured.



Standard tip

12,5 cm
Inox
0,4 x 0,2 mm

No. 11210-20

Standard tip

12,5 cm
Epoxy coated inox
0,4 x 0,2 mm

No. 11210-10

Standard tip

13,5 cm
Inox
0,20 x 0,12 mm

No. 11203-23

Angled 45°

13,5 cm
Inox
0,20 x 0,12 mm

No. 11203-25

Biology tip

13,5 cm
Dumoxel
0,10 x 0,06 mm

No. 11200-33

FORCEPS

For information on
Dumont alloys
and tip styles,
see page 209.

Dumont Mini Forceps

FORCEPS

For information on
Dumont alloys
and tip styles,
see page 209.



Style 3

7 cm

Inox

0,12 x 0,07 mm

No. 11200-10



Style 4

7 cm

Inox

0,14 x 0,06 mm

No. 11200-12



Style 5


8 cm

Inox

0,09 x 0,05 mm

No. 11200-14

Dumont #5AC



300%

Biology Tip
Type #5 Anti-Capillary

11,5 cm

Inox

0,07 x 0,02 mm

No. 11258-20

O-Rings for Forceps

These rubber O-rings fit our forceps
and keep them closed with variable
degrees of pressure on the tips.

NOT
AUTOCLAVABLE

Pkg. of 10

No. 11200-00



Gold Plated Forceps

Gold-Plating gives corrosion resistance to most chemicals, salts and acids. These forceps are perfect for any application where corrosion resistance and toughness are primary requirements.



Straight

11,5 cm
0,1 x 0,06 mm

Gold Plated
Stainless Steel

No. 11915-11



Curved

11,5 cm
0,2 x 0,1 mm

Gold Plated
Stainless Steel

No. 11917-12



Angled 30°

11,5 cm
0,1 x 0,06 mm

Gold Plated
Stainless Steel

No. 11911-11

Diamond Coated Forceps

Diamond coating provides increased corrosion resistance and greater protection from wear and tear, thus extending the operational life of the forceps, as well as providing a stronger grip at the tips.



Straight

14 cm
0,15 x 0,1 mm

Diamond Coated
Stainless Steel

No. 11990-14



Angled 30°

14 cm
0,15 x 0,1 mm

Diamond Coated
Stainless Steel

No. 11991-14



Straight

12 cm
1,65 x 0,05 mm

Diamond Coated
Stainless Steel

No. 11992-12

Mirror Finish Broad Shanks

These forceps have long and narrow shanks with thin, highly polished tips to prevent tissue from sticking. Ideal for embryos, nerves, or other adhesive tissues.



Straight

13 cm
0,3 x 0,15 mm

No. 11923-13

Forceps with Replaceable Tips

Stainless steel body with replaceable plastic tips. Tips are made of polyvinylidene fluoride (PVDF) and are electrostatic discharge (ESD) safe. Body and plastic tips are corrosion and abrasion resistant.



Straight

13 cm
0,5 x 0,5 mm

Stainless Steel body with
replaceable plastic tips

No. 11980-13



Pkg. of 2 tips and 3 screws
Replacement plastic tips

Straight, strong

4 cm
0,5 mm

No. 11988-01

Pkg. of 2 tips and 3 screws
Replacement plastic tips

Straight, very fine

4 cm
0,5 mm

No. 11988-02

Pkg. of 2 tips and 3 screws
Replacement plastic tips

Curved

4 cm
0,5 mm

No. 11988-03

Fine Forceps Broad Shanks

Straight forceps with sturdy coarse tips.



12 cm
0,2 x 0,12 mm
No. 11445-12

Fine Forceps Mirror Finish

The last 5 mm on the inside of the tips are highly polished to help prevent tissue from sticking. Ideal for embryos, nerves, or other adhesive tissue.



11 cm
0,10 x 0,06 mm
No. 11412-11



11 cm
0,05 x 0,02 mm
No. 11413-11

Fine Forceps Self-Closing

Ideal for lifting and holding EM grids. Consistent pressure at tips. Correctly tensioned "cross-over" shanks come to a microscopically matched point.



11,5 cm
0,08 x 0,04 mm
No. 11480-11



11 cm
0,10 x 0,06 mm
No. 11485-11



11,5 cm
0,20 x 0,12 mm
No. 11487-11

Moria Ultra Fine Forceps

These forceps have the most finely executed points available. They are extremely light, superbly balanced, and have a gentle spring action. The shanks are specially designed with an arch which ensures maximum force to the tips. The tips close flat for 6 mm from the points to give a large holding area. The tips have been expertly finished under magnification.



Moria MC40



Straight

11 cm

No. 11370-40

Moria MC40/B



Curved

11 cm

No. 11370-42

Moria 9980



Straight

13 cm

No. 11399-80

Moria 9987



Curved

13 cm

No. 11399-87

Moria Iris Forceps

Balanced with gentle spring action and tips that are perfectly matched.



Moria MC32/B

Smooth
Blunt10 cm
Tip width: 0,5 mm

No. 11373-22

Moria MC31/B

Smooth
Blunt10 cm
Tip width: 0,5 mm

No. 11373-12

Moria MC32



Microserrated

No. 11370-32

Moria MC31



Microserrated

No. 11370-31

Fine Forceps

Straight

11 cm
0,3 x 0,25 mm

No. 00108-11




Angled 45°

11 cm
0,3 x 0,25 mm

No. 00109-11

SuperGrip Forceps

The tips are processed with MicroCrystals®, a method developed by S&T® to infuse tungsten carbide on stainless steel. The formation of this ultra-fine coating guarantees atraumatic handling of tissue while improving the grip dramatically. The recommended use is where forceps are utilized as needle holders or when slippery tissues might be encountered.



300%

The gripping surface is similar to a match head.



Straight

11 cm
0,3 x 0,25 mm


No. 00632-11




Angled 45°

11 cm
0,3 x 0,25 mm

No. 00649-11

Bonn

20x magnification




Smooth

7 cm

Tip width: 0,3 mm

Tying Platform: 5 mm

No. 11083-07




1x2 teeth

7 cm


Tip width: 0,45 mm

Teeth width: 0,1 mm

No. 11084-07

Mirror Finish Bonn Forceps

20x magnification



Straight

1x2 teeth

10 cm

Tip width: 0,12 mm


Tying Platform: 5 mm

Partly polished


No. 11090-10

Mirror Finish Suture Tying Forceps with Teeth

A unique tip style fusing a suture tying platform and a 1x2 teeth tip for tissue manipulation.



20x magnification



Straight


9 cm

Tip width: 0,3 mm

Platform with 1x2 teeth

No. 11090-09


Mirror Finish Suture Tying Forceps



Straight
10 cm
Tip diameter: 0,3 mm
Tying Platform: 5 mm
Partly polished
No. 11093-10

McPherson Forceps

MOST
POPULAR




Angled 45°
Smooth
9 cm
Tip width: 0,4 mm
Tying Platform: 6 mm
No. 11063-07

Suture Tying Forceps

The platform tips close flat. This ensures a reliable grip on sutures.


S&T®



13,5 cm
Tip diameter: 0,3 mm
Tying platform: 6 mm
No. 00272-13

Suture Tying Forceps


The platform tips close flat. This ensures a reliable grip on sutures.



10 cm
Tip diameter: 0,5 mm
Tying platform: 5 mm
No. 18025-10

**Round Handled
Suture Tying Forceps**



Round handles allow easy fingertip adjustments.



10,5 cm
Tip diameter: 0,35 mm
Tying platform: 5,5 mm
No. 18026-10

Vial Forceps

Tips taper for 6 cm allowing easy access to specimens in vials.



15 cm
Tip diameter: 0,5 mm
No. 11006-15

Vitro Grasping Forceps


14 cm
Serrated
Tip diameter: 0,7 mm
No. 11900-14

Graefe ForcepsStraight
Serrated


7 cm
Tip width: 0,6 mm
No. 11064-07

Curved
Serrated

7 cm
Tip width: 0,6 mm
No. 11065-07

Straight
1x2 teeth


No. 11066-07

Curved
1x2 teeth

No. 11067-07

Bishop-Harman

Extra light, fenestrated handles ensure a good grip.




10x magnification



Serrated

8,5 cm
Tip width: 0,5 mm
No. 11069-08



1x2 teeth

8,5 cm
Tip width: 0,5 mm
No. 11070-08

Octagon Forceps


The octagon shaped grip provides a particularly secure finger hold.

Straight
Serrated


9,5 cm
Tip width: 0,8 mm
No. 11040-08

Curved
Serrated

9,5 cm
Tip width: 0,8 mm
No. 11041-08

Straight
1x2 teeth

No. 11042-08

Curved
1x2 teeth

No. 11043-08

Extra Fine Graefe Forceps

80% of
actual size

Serrated

10 cm

Tip width: 0,5 mm

No. 11150-10



Serrated

10 cm

Tip width: 0,5 mm

No. 11151-10



Serrated

10 cm

Tip width: 0,5 mm

No. 11152-10



Smooth

No. 11158-10




Smooth

No. 11159-10



Smooth

No. 11160-10



1x2 teeth

No. 11153-10



1x2 teeth

No. 11154-10



1x2 teeth

No. 11155-10

Graefe Forceps

MOST
POPULAR11051-10
11054-10

•TITANIUM•

11650-10 – 11652-10
For more information
see page 203.

Serrated

Tip width: 0,8 mm

No. 11050-10



Serrated

Tip width: 0,8 mm

No. 11051-10



Serrated

Tip width: 0,8 mm

No. 11052-10



Serrated

Tip width: 0,8 mm

No. 11049-10

Titanium

Tip width: 1,2 mm

No. 11650-10

Titanium


Tip width: 1,2 mm

No. 11651-10

Titanium

Tip width: 1,2 mm

No. 11652-10



1x2 teeth

Tip width: 0,8 mm

No. 11053-10



1x2 teeth

Tip width: 0,8 mm

No. 11054-10



1x2 teeth

Tip width: 0,8 mm

No. 11055-10

Sliding Lock Forceps

Smooth

11 cm

Tip width: 1 mm
Slide lock

No. 11062-10

Episcleral Forceps

The tips of these forceps arch toward each other, closing flat. Exceptional gripping power is possible with minimal damage to tissue. A limit pin prevents excessive applied pressure.



300%



10 cm

Tip thickness: 0,2 mm
Tip Dimension: 1,3 x 3 mm

No. 11071-10

Pierse Fixation Forceps

For exceptional gripping power with no damage to sclera, vascular sheath and epineurium.



1800%

Tissue retention
hole of 0,3 mm
behind tips.

Straight

13 cm

Tip Dimension: 0,1 x 0,4 mm

No. 18155-13

Plastic Forceps


These forceps are non-magnetic, acid resistant and able to withstand temperatures up to 240°C (464°F). Tips do not spread under pressure. Plastic forceps are made of Acetyl and will be sold in packages of 10 pieces.



Sharp
Straight
11,5 cm
0,40 x 0,25 mm
No. 11700-00




Blunt
Straight
11,5 cm
1,5 x 0,40 mm
No. 11700-03



Sharp
Curved
11,5 cm
0,30 x 0,20 mm
No. 11700-01




Blunt
Curved
11,5 cm
1,5 x 0,40 mm
No. 11700-04



Sharp
Angled 45°
11,5 cm
0,30 x 0,20 mm
No. 11700-02


Cross Action Forceps

These tissue transport forceps have a gentle action.



Cross serrated
10 cm
Tip Dimension: 4 x 8 mm
No. 11075-00

Cover Glass Forceps



Straight	Angled
10,5 cm	10,5 cm
Tip Dimension: 4 x 5 x 10 mm	Tip Dimension: 4 x 5 x 10 mm
No. 11074-02	No. 11073-10

Ring Forceps

For transport of tissue and organs. With limit screw to prevent excessive applied pressure.



8,5 cm
Inner Diameter: 1 mm
Outer Diameter: 1,5 mm

No. 11101-09



9 cm
Inner Diameter: 2,2 mm
Outer Diameter: 3 mm

No. 11103-09



9 cm
Inner Diameter: 4,8 mm
Outer Diameter: 6 mm

No. 11106-09

Circle Forceps

Cupped forceps for atraumatic gripping of organs and tissues.



11,5 cm
Tip diameter: 0,7 mm

No. 11072-11

Live Insect Handling Forceps

Very delicate.



Narrow
10 cm

No. 26029-10




Wide
10 cm

No. 26030-10

Micro-Adson

Our most delicate version of the popular Adson forceps.



Serrated

12 cm
Tip width: 0,8 mm

No. 11018-12


1x2 teeth

12 cm
0,8 mm

No. 11019-12

Micro-Adson with Fenestrated Handle

Our delicate Micro-Adson forceps, with five windows to aid in stability and easy finger adjustments.



Serrated

12 cm
0,8 mm

No. 11020-12

1x2 teeth

12 cm
0,5 mm

No. 11032-12

1x2 teeth with
6 mm platform

12 cm
0,8 mm

No. 11022-12

Adson


Wide, strong grip with narrow tips.

**MOST
POPULAR**

11006-12
11027-12

**TUNGSTEN
CARBIDE**

11506-12
For more information
see page 203.



Serrated

12 cm
1,5 mm

No. 11006-12

1x2 teeth



12 cm
1,5 mm

No. 11027-12

Cross serrated
tungsten carbide
inlaid tip

12 cm
1,5 mm

No. 11506-12

Adson DeBakey**Atraumatic
DeBakey Serrations**

This tooth pattern, as shown, consists of one row of atraumatic teeth meshed with two opposite rows of atraumatic teeth.


These two forceps are ideal for gentle, atraumatic gripping of vessels, tissue and other delicate structures.



DeBakey serrations

12,5 cm
Tip width: 1,5 mm

No. 11617-12



DeBakey serrations

15 cm
Tip width: 1,5 mm


No. 11616-15

Adson-Graefe

Fine sharp teeth

12,5 cm
Tip width: 4,5 mm

No. 11030-12

Adson-Brown

Shark teeth

12 cm
Tip width: 1,4 mm

No. 11627-12


Standard Pattern

Broad pattern with strong spring action.



Straight
Serrated

Tip width (mm)	Length (cm)	Item
2,5	12	No. 11000-12
2,8	13	No. 11000-13
3,25	14,5	No. 11000-14
3,3	16	No. 11000-16
3,7	18	No. 11000-18
3,75	20	No. 11000-20
3,9	25	No. 11000-25
4	30	No. 11000-30



Curved
Serrated

Tip width (mm)	Length (cm)	Item
2,5	12	No. 11001-12
2,8	13	No. 11001-13
3,25	14,5	No. 11001-14
3,3	16	No. 11001-16
3,7	18	No. 11001-18
3,75	20	No. 11001-20

Narrow Pattern

Lighter spring action than standard pattern
for grasping more delicate tissue.


•TITANIUM•

11602-14
For more information
see page 203.



Straight
Serrated

Tip width (mm)	Length (cm)	Item
2	12	No. 11002-12
2,15	13	No. 11002-13
2,25	14,5	No. 11002-14
2,25	14,5	No. 11602-14
2,35	16	No. 11002-16
2,4	18	No. 11002-18
2,75	20	No. 11002-20
3,2	25	No. 11002-25



Curved
Serrated

Tip width (mm)	Length (cm)	Item
2	12	No. 11003-12
2,15	13	No. 11003-13
2,25	14,5	No. 11003-14
2,35	16	No. 11003-16
2,4	18	No. 11003-18
2,75	20	No. 11003-20

Semken

Slim pattern with long, narrow shanks and gentle spring action.

•TITANIUM•

11608-15

For more information
see page 203.



Straight
Serrated

Tip width (mm)	Length (cm)	Item
1,3	13	No. 11008-13
1,5	15	No. 11008-15
1,5	15	No. 11608-15



Curved
Serrated

Tip width (mm)	Length (cm)	Item
1,3	13	No. 11009-13
1,5	15	No. 11009-15

Mouse Handling Forceps

Long handles make picking up small rodents easy and help prevent bites.



21 cm
Forceps with serrated silicone tips
Gripping surface: 15 x 20 mm

No. 11036-20



Replacement silicone tips
Pkg of 2 tips

No. 11036-90

Tissue Forceps



Straight
Standard 1x2 teeth
12 cm

Straight Narrow 1x2 teeth 10 cm No. 11023-10	Straight Narrow 1x2 teeth 14,5 cm No. 11023-14
12 cm No. 11023-12	15,5 cm No. 11023-15


No. 11021-12 14,5 cm No. 11021-14 15,5 cm No. 11021-15	Straight Standard 2x3 teeth 14,5 cm No. 11022-14 15,5 cm No. 11022-15
--	--

Tissue Forceps

Used with rongeur No. 16021-14
(page 36) in laminectomy to securely
hold spinal column.



Curved
Narrow 1x2 teeth
14 cm
No. 11029-14

Gillies


Cross-Serrated
1x2 teeth

15 cm
Tip Dimension: 1,5 x 2 mm

No. 11028-15

Russian

Special cupped pattern
for superior grip.




15 cm
Tip Dimension: 6 x 11 mm

No. 11026-15

Organ Holding Forceps


With limit screw to prevent
excessive applied pressure.



Serrated
14,5 cm
Tip Dimension: 8 x 11 mm

No. 11031-15

London Forceps



Angled
Serrated

16 cm

No. 11080-02

Titanium Forceps

•TITANIUM•

For more information
see page 203.




Straight
Serrated

16 cm

No. 11602-16


Gerald Forceps



Straight
Serrated

18 cm

No. 11014-18





Curved
Serrated

18 cm

No. 11015-18

Taylor Forceps

Serrated and smooth versions available.

**Potts-Smith Forceps**

Ultra Fine Hemostat



Straight
Smooth

12 cm
Tip width: 0,6 mm

No. 13020-12



Curved
Smooth

12 cm
Tip width: 0,6 mm

No. 13021-12

Fine Hemostat



Straight
Smooth

12 cm
Tip width: 1 mm

No. 13006-12



Curved
Smooth

12 cm
Tip width: 1 mm

No. 13007-12

Micro-Mosquito



Straight
Serrated

12 cm
Tip width: 1,3 mm


No. 13010-12



Curved
Serrated

12 cm
Tip width: 1,3 mm

No. 13011-12



Angled 90°
Serrated

12 cm
Tip width: 1,3 mm

No. 13012-12

Halsted-Mosquito



13108-12



13008-12



13109-12



Straight
Serrated

12,5 cm
Tip width: 1 mm

No. 13008-12



Curved
Serrated

12,5 cm
Tip width: 1 mm

No. 13009-12



Straight
Serrated

12,5 cm
Tip width: 1 mm
Left handed

No. 13108-12




Curved
Serrated

12,5 cm
Tip width: 1 mm
Left handed

No. 13109-12


Crile



Straight
Serrated

14 cm
Tip width: 1,5 mm

No. 13004-14



Curved
Serrated

14 cm
Tip width: 1,5 mm

No. 13005-14


Kelly



Straight
Serrated

14 cm
Tip width: 1,5 mm

No. 13018-14



Curved
Serrated


14 cm
Tip width: 1,5 mm

No. 13019-14

Baby-Mixter

Providence-Hospital

HEMOSTATS



Curved
Serrated
14 cm
Tip width: 1,3 mm
No. 13013-14



Straight
Serrated
14 cm
Tip width: 1,5 mm
No. 13014-14




Curved
Serrated
14 cm
Tip width: 1,5 mm
No. 13015-14

Hartman

**MOST
POPULAR**

13003-10



**Straight
Serrated**

10 cm
Tip width: 1 mm

No. 13002-10



**Curved
Serrated**

10 cm
Tip width: 1 mm

No. 13003-10

Backhaus Towel Clamp



9 cm

No. 11095-09



13 cm

No. 11095-13

Allis Clamp

Tissue and organ holding forceps.



Allis



5x6 teeth
15 cm
Tip width: 4 mm
No. 11091-15

Allis-Baby



4x5 teeth
12 cm
Tip width: 3 mm
No. 11092-12

Lorna Clamp



Serrated
10 cm
Tip width: 6 mm
No. 11096-10



Serrated
13 cm
Tip width: 7 mm
No. 11096-13