

**Pin Holder**

Light weight



90% of actual size

17 cm  
Nickel-plated chuck  
Aluminum handle  
Jaw opening: 0–1 mm

No. 26018-17

**Pin Holder**




Moria MC1

12 cm  
Nickel-plated  
Jaw opening: 0–1 mm

No. 26016-12

**Self Holding Pin and Bone Screw Forceps**

Designed with a serrated tip and grooved cup for a stable and secure grip, these forceps hold pins and screws allowing precise placement at your target area. They are compatible with all Fine Science Tools insect pins, minuten pins and bone screws.



200%





14 cm  
Polished stainless steel

No. 26100-00

**Insert Pins**

All insert pins fit any of our pin holders.  
Pins are antimagnetic and acid resistant.



**Black Anodized Minutien Pins**

These high quality pins are black for best visibility under high magnification and bright light.

Tip diameter: 0,0175 mm  
Rod diameter: 0,15 mm  
10 mm

Box of 500 pins  
No. 26002-15 **MOST POPULAR**



**Stainless Steel Minutien Pins**

The thinnest and finest pins available.

Tip diameter: 0,0125 mm  
Rod diameter: 0,1 mm  
10 mm


Box of 500 pins  
No. 26002-10 **MOST POPULAR**

Tip diameter: 0,02 mm  
Rod diameter: 0,2 mm  
10 mm

Box of 500 pins  
No. 26002-20 **MOST POPULAR**

**Black Anodized Insect Pins**

Nylon heads  
Tip diameter: 0,03 mm  
Pkg. of 100 pins



Size	000	00	0	1	2	3	4	5	6	7
Rod diameter (mm)	0,25	0,30	0,35	0,40	0,45	0,50	0,55	0,60	0,65	0,70
No. 26000	-25	-30	-35	-40	-45	-50	-55	-60	-65	-70

**Stainless Steel Insect Pins**

Nylon heads  
Tip diameter: 0,03 mm  
Pkg. of 100 pins

Size	000	00	0	1	2	3	4	5	6	7
Rod diameter (mm)	0,25	0,30	0,35	0,40	0,45	0,50	0,55	0,60	0,65	0,70
No. 26001	-25	-30	-35	-40	-45	-50	-55	-60	-65	-70

**Etched Tungsten Micro Needles**

These tungsten micro needles are drawn to a 1 micron tip and three different rod diameters are offered. Each provides lateral stability as well as minimal spring action during manipulation. The micro needles can be used with our pin holders shown on page 112.



Tip diameter: 1 µm (0,001 mm)  
12 mm

Box of 10 needles

Rod diameter (mm)	Item
0,125	No. 10130-05
0,25	No. 10130-10
0,5	No. 10130-20