

**PHYMEP**

21 RUE CAMPO-FORMIO 75013 PARIS - tél : 01 44 24 25 18 - fax : 01 44 24 23 57

email : phymep@wanadoo.frsite : <http://www.phymep-sarl.com/>**A-M Systems****FIL NU**

Catalog Number	Diameter		Nominal AWG	Length / Spool	
	in	µm		ft	m

GOLD / OR

750020	0.002	50.800	44	10	3.048
750000	0.003	76.200	41	10	3.048

SILVER / ARGENT

781000	0.003	76.200	41	25.000	7.620
781500	0.005	127	36	25.000	7.620
782000	0.008	203.200	32	25.000	7.620
782500	0.010	254.000	30	25.000	7.620
783000	0.015	381.000	27	25.000	7.620
783500	0.025	635.000	22	25.000	7.620

PLATINIUM / PLATINE

764500	0.001	25.400	50	10	3.048
765000	.002	50.800	44	10	3.048
766000	0.003	76.200	41	10	3.048
767000	0.005	127	36	10	3.048

PLATINIUM IRIDIUM / PLATINE IRIDIÉ

767400	0.001	25.400	50	10	3.048
767500	0.002	50.800	44	10	3.048
767600	0.003	76.200	41	10	3.048
767700	0.005	127	36	10	3.048

TUNGSTEN / TUNGSTENE

797500	0.002	50.800	44	25.000	7.620
797550	0.003	76.200	40	25.000	7.620
797600	0.005	127.000	36	25.000	7.620
797650	0.008	203.200	32	25.000	7.620



STAINLESS STEEL / ACIER INOX

Catalog Number	Diameter		Nominal AWG	Hardness		Length / Spool	
	in	µm		Strength	kPSI	ft	m
793600	0.002	50.800	44	Annealed	90-120	100.000	30.480
792600	0.002	50.800	44	Half	155-185	100.000	30.480
794600	0.002	50.800	44	Full	215-245	100.000	30.480
793700	0.003	76.200	41	Annealed	90-120	100.000	30.480
792700	0.003	76.200	41	Half	155-185	100.000	30.480
794700	0.003	76.200	41	Full	215-245	100.000	30.480
792800	0.005	127.000	38	Annealed	90-120	100.000	30.480
793800	0.005	127.000	38	Half	155-185	100.000	30.480
794800	0.005	127	38	Full	215-245	100.000	30.480
792900	0.008	203.200	32	Annealed	90-120	100.000	30.480
793900	0.008	203.200	32	Half	155-185	100.000	30.480
794900	0.008	203.200	32	Full	215-245	100.000	30.480
793000	0.010	254.000	30	Annealed	90-120	100.000	30.480
794000	0.010	254.000	30	Half	155-185	100.000	30.480
795000	0.010	254.000	30	Full	215-245	100.000	30.480
794100	0.020	508.000	24	Annealed	90-120	100.000	30.480
793100	0.020	508.000	24	Half	155-185	100.000	30.480
795100	0.020	508.000	24	Full	215-245	100.000	30.480