








Ear Bars

Set #	Model #	Name	Description	Tip Shape Shown
1	821 921 1721	Non-Rupture (Mouse, Neonatal Rat)	Zygoma ear cups (serrated and soft tissue)	
2	822 922 1722	Non-Rupture (Mouse, Siberian Hamster)	60° Non-Rupture tip <i>Hollow Ear Bars are also available</i>	
3	851 951 1751	Standard Solid (Cat / Monkey)	20° tapered tip with a 1.5 mm radius	
4	852 952 1752	Non-Rupture (Cat / Monkey)	20° tapered tip to a 4.8 mm shoulder with a 3 mm dia. by 2 mm long protrusion	
5	853 953 1753	Hollow (Cat / Monkey)	4 mm hole through for auditory stimulation	
6	855 955 1755	Non-Rupture (Bird, Guinea Pig, Prairie Vole, Rat, Squirrel Monkey)	45° Non-Rupture tip with a 0.8 mm radius <i>Hollow Ear Bars are also available</i>	
7	856 956 1756	(Pigeon)	3 mm dia. by 6 mm long tip with a full radius on end	
8	857 957 1757	(Rat)	18° tip with 0.8 mm radius <i>Hollow Ear Bars are also available</i>	