



▶ **Micromanipulator MD 4** with a second tool holder

Specifications

Travel range:	X-axis: 37 mm Y-axis: 20 mm Z-axis: 25 mm fine adjustment: 10 mm
Weight:	700 g
Packaging size / -weight:	40 x 32 x 21 cm / 2 kg

Order Information

MD 4 right with clamp 1/2"	Part No.: 00-42-105-0000
MD 4 left with clamp 1/2"	Part No.: 00-42-106-0000
MD 4 right with tilting device	Part No.: 00-42-121-0000
MD 4 left with tilting device	Part No.: 00-42-122-0000

Accessories

Additional Z-axis stop	Part No.: 00-45-201-6000
Magnetic stand M1	Part No.: 00-49-901-6000

Features

The MD 4 is a compact micromanipulator for manual movement in three axes (X, Y, and Z). It is equipped with a second tool holder. This additional tool holder cannot only be positioned in X- and Y-direction, independent from the micromanipulator, but may also be tilted and swiveled via two fine adjustment screws.

The scales mounted on its guide ways allow read off of the coarse adjustment with an accuracy of 0,1 mm. The additional fine adjustment is operated by turning the micrometer screw, the read off accuracy is 1/100 mm.

Since all control knobs are mounted in line on the rear (not on the side), a number of user friendly MD 4 micromanipulators can be set up right next to one another.

For the purpose of avoiding that the Z-axis will descend when working with weighty loads (drift), there is an optional break (Z-axis stop) available.

Märzhäuser Wetzlar GmbH & Co. KG
In der Murch 15
D-35579 Wetzlar

Tel.: +49 (0) 64 41/91 16-0
Fax: +49 (0) 64 41/91 16-40
eMail: info@marzhauser.com