





MMJ right-sided
with clamp 1/2"




MMJ left-sided
with clamp 1/2"

Manual Micromanipulator MMJ with Mechanical Joystick

Specifications

Travel range	X-axis: 37 mm Y-axis: 20 mm Z-axis: 25 mm
Travel range joystick	0,35 to 4,0 mm
Material	aluminium
Surface	anodic coating, black lacquered
Weight	950 g



Order Information

Order No.: 00-42-107-0000 (right-sided, with clamp 1/2")
Order No.: 00-42-108-0000 (left-sided, with clamp 1/2")

Product Features

- ▶ manual movement via control knobs in 3 axes (XYZ)
- ▶ additional mechanical joystick for XY with travel range from 0,35 to 4,0 mm at full detection
- ▶ coarse and fine adjustment of Z-axis, ratio 1 : 10
- ▶ tool holder for holding a piezo translator (e. g. PM 20)
- ▶ tool holder tiltable up to 90°
- ▶ stop screw assures that tool returns to focus level after descending

Accessories

Magnetic stand M1
Order No.: 00-49-901-6000

Magnetic stand M2 (rotatable/tiltable)
Order No.: 00-49-902-6000