



PHYMEP

21 RUE CAMPO-FORMIO 75013 PARIS - tél : 01 44 24 25 18 - fax : 01 44 24 23 57

email : phymep@wanadoo.fr

site : <http://www.phymep-sarl.com/>

Vis Inox



Diamètre Tête 2.50 mm

Diamètre Tige 1.57 mm

RAT – Plastics One

Réf.	Longueur Tige	Diamètre mèche	Tournevis
0-80 x 1/16	1.6 mm	D#56 – 1.18 mm	SD-80
0-80 x 1/8	3.2 mm	D#56 – 1.18 mm	SD-80
0-80 x 3/16	4.8 mm	D#56 – 1.18 mm	SD-80
0-80 x 3/32	2.4 mm	D#56 – 1.18 mm	SD-80

SOURIS – Plastics One

Réf.	Longueur Tige	Diamètre mèche	Tournevis
00-96x1/16	1.6mm	D#60 – 1.0 mm	SD-96
00-96 x 3/32	2.4mm	D#60 – 1.0 mm	SD-96

Vis Auto Perforante

Réf.	Longueur Tige	Diamètre Tige
19010-00	4.7mm	1.17mm
19010-10	4mm	0.85mm



Vis Tête Plate

