

OD mm	ID mm	Wall Thickness	Tolerance mm	Length	Qty/Pkg	Model	Order #
1.0	0.50	0.5	± 0.05	75 mm	500	GC100FS-7.5	30-0034
				100 mm	500	GC100FS-10	30-0032
				150 mm	500	GC100FS-15	30-0033
	0.58	0.42		75 mm	500	GC100F-7.5	30-0021
				100 mm	500	GC100F-10	30-0019
				150 mm	500	GC100F-15	30-0020
1.2	0.69	0.51	± 0.05	75 mm	350	GC120F-7.5	30-0046
				100 mm	350	GC120F-10	30-0044
				150 mm	350	GC120F-15	30-0045
1.5	0.86	0.64	± 0.05	75 mm	225	GC150F-7.5	30-0060
				100 mm	225	GC150F-10	30-0057
				150 mm	225	GC150F-15	30-0058
2.0	1.16	0.42	± 0.1	75 mm	125	GC200F-7.5	30-0076
				100 mm	125	GC200F-10	30-0074
				150 mm	125	GC200F-15	30-0075
3.0	1.62	1.38	± 0.1	75 mm	55	GC300F-7.5	30-0084
				100 mm	55	GC300F-10	30-0082
				150 mm	55	GC300F-15	30-0083