

OD mm	ID mm	Wall Thickness mm	Tolerance mm	Length	Qty/Pkg	Model	Order #
1.0	0.58	0.21	± 0.05	75 mm	500	GC100-7.5	30-0018
				100 mm	500	GC100-10	30-0016
				150 mm	500	GC100-15	30-0017
1.2	0.69	0.255	± 0.05	75 mm	350	GC120-7.5	30-0043
				100 mm	350	GC120-10	30-0042
				150 mm	350	GC120-15	30-0041
1.5	0.86	0.32	± 0.05	75 mm	225	GC150-7.5	30-0056
				110 mm	225	GC150-11	30-0054
				100 mm	225	GC150-10	30-0053
				150 mm	125	GC150-15	30-0055
2.0	1.16	0.42	± 0.1	75 mm	125	GC200-7.5	30-0073
				100 mm	125	GC200-10	30-0070
				150 mm	125	GC200-15	30-0071
3.0	1.62	1.38	± 0.1	75 mm	55	GC300-7.5	30-0127
				100 mm	55	GC300-10	30-0080
				150 mm	55	GC300-15	30-0081