

Ultra Fine Hemostat

These are the finest hemostats in our product line.



Straight Smooth

12 cm
Tip width: 0.6 mm

No. 13020-12



Curved Smooth

12 cm
Tip width: 0.6 mm

No. 13021-12

Fine Hemostat



Straight Smooth

12 cm
Tip width: 1 mm

No. 13006-12



Curved Smooth

12 cm
Tip width: 1 mm

No. 13007-12

Micro-Mosquito

Mosquito



13108-12
13109-12



HEMOSTATS



Straight Serrated

12.5 cm
Tip width: 1 mm


No. 13010-12



Curved Serrated

12.5 cm
Tip width: 1 mm

No. 13011-12



Angled 90° Serrated

12.5 cm
Tip width: 1 mm

No. 13012-12



Straight Serrated

12.5 cm
Tip width: 1.1 mm

No. 13008-12

Left handed

No. 13108-12



Curved Serrated

12.5 cm
Tip width: 1.1 mm


No. 13009-12

Left handed

No. 13109-12

Crile

Kelly



Straight Serrated

14 cm
Tip width: 1.5 mm

No. 13004-14



Curved Serrated

14 cm
Tip width: 1.5 mm

No. 13005-14



Straight Serrated

14 cm
Tip width: 1.5 mm

No. 13018-14




Curved Serrated

14 cm
Tip width: 1.5 mm

No. 13019-14

Adson-Baby

Providence-Hospital



**Curved
Serrated**

14 cm
Tip width: 1.3 mm

No. 13013-14



**Straight
Serrated**

14.5 cm
Tip width: 1.5 mm

No. 13014-14




**Curved
Serrated**

14.5 cm
Tip width: 1.5 mm

No. 13015-14

HEMOSTATS

Hartmann



Straight Serrated

10 cm
Tip width: 1 mm

No. 13002-10



Curved Serrated

10 cm
Tip width: 1 mm

No. 13003-10

Backhaus Towel Clamp



Curved Round tips

9 cm

No. 11095-09




Curved Round tips

13 cm

No. 11095-13

Allis


Tissue and organ holding forceps.



**Straight
5x6 Teeth**

15 cm

No. 11091-15



**Straight
4x5 Teeth**

12 cm

No. 11092-12


Lorna Towel Clamp



**Curved
Serrated**

10 cm
Tip width: 7 mm

No. 11096-10



**Curved
Serrated**

13 cm
Tip width: 7 mm

No. 11096-13